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

2. Logga in på myAbilia

The desktop login form features the myABILIA logo at the top. Below it, there are two input fields: one for 'Användarnamn/E-post' (Username/Email) and one for 'Lösenord' (Password). A red 'Logga in' button is positioned below the password field. At the bottom, there is a link that says 'Glömt ditt lösenord?' (Forgot your password?).

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


- Bild:** A dashed box with a plus icon and a callout '1' for selecting an image.
- Namn:** A text input field with a callout '2' for entering the activity name.
- Datum:** A date picker showing '2022-11-02' with a callout '3' for selecting the date.
- Starttid:** A time picker showing '10' with a callout '4' for selecting the start time.
- Sluttid:** A time picker showing '10' with a callout '4' for selecting the end time.
- Heldag:** A checkbox with a callout '4' for selecting 'Heldag'.
- Förinställningar:** A section with three options: 'Ta bort efteråt' (callout '5'), 'Kvittera' (callout '6'), and 'Höger' (callout '7').
- Synlig för:** A section with three options: 'Bara mig' (callout '8'), 'Alla mina stödpersoner' (callout '8'), and 'Utvalda stödpersoner' (callout '8').


At the bottom left, there is a button 'Avbryt' (Cancel).

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The screenshot shows the 'MEMOplanner' app interface for setting a new activity. The title is 'Ny aktivitet' and there is a 'Basaktiviteter' button in the top right. The interface is divided into four tabs: 'Namn & Bild', 'Larm & Påminnelser', 'Upprepning', and 'Kopplade funktioner'. The 'Larm & Påminnelser' tab is active and contains three sections:

- Larm (1):** Options for alarm type: 'Larm + vibration' (selected), 'Endast larm', 'Endast vibration', 'Tyst larm', and 'Inget larm'. There is also a checkbox for 'Larm endast vid starttid'.
- Påminnelse (2):** Time intervals for reminders: '5 minuter', '15 minuter', '30 minuter', '1 timme', and '2 timmar'.
- Tal (3):** Options for voice recording: 'Tal vid starttid', 'Tal vid sluttid', and a 'Starta inspelning' button. A small inset window shows the 'Tal' settings with 'Tal vid sluttid' checked and a recording progress bar from 00:00 to 03:00.

An 'Avbryt' button is located at the bottom center of the screen.

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
The screenshot displays the 'Stödpersoner' (Support persons) section of the ABILIA system. A sidebar on the left contains navigation options, with 'Stödpersoner' highlighted. The main content area shows a table with columns for 'Namn', 'Roll', 'E-postadress', and 'Hantera'. A modal window titled 'Bjud in stödperson' is open, showing a form for adding a new support person. The form includes a field for 'E-postadress' with the value 'epost.stodperson@domain.se'. A green 'Bjud in' button is located at the bottom right of the modal, and an 'Avbryt' (Cancel) button is at the bottom left.

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