


PHILIPS
RESPIRONICS

InnoSpire Deluxe

InnoSpire Deluxe

Index

English	EN
Български	BG
Česky	CZ
Hrvatski	HR
Dansk	DA
Deutsch	DE
Español	ES
Eesti	EE
Français	FR
Ελληνικά	EL
Indonesia	ID
Italiano	IT
Latviski	LV
Lietuviškai	LT
Magyar	HU
Nederlands	NL
Norsk	NO
Polski	PL
Português	PT
Română	RO
Srpski	SR
Slovenčina	SK
Slovenščina	SL
Suomi	FI
Svenska	SV
ภาษาไทย	TH
Tiếng Việt	VI
Türkçe	TR
العربية	AR
한국어	KO
עברית	IL
EMC data sheets (EN)	

4
10
16
22
28
34
40
46
52
58
64
70
76
82
88
94
100
106
112
118
124
130
136
142
148
154
160
166
172
178
184
190

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A. Instructions for use

- Innospire Deluxe needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in this document;
- Portable and mobile RF communications equipment can affect Innospire Deluxe.

B. Technical description

- Warning** that the use of **accessories, transducers and cables** other than those specified with the exception of transducers and cables sold by the manufacturer of the Innospire Deluxe as replacement parts for internal components, may result in increased EMISSIONS or decreased IMMUNITY of the Innospire Deluxe.
- Warning** that the Innospire Deluxe should not be used adjacent to or stacked with other equipment.

Guidance and manufacturer's declaration – electromagnetic emissions		
Emissions	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The Innospire Deluxe uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Innospire Deluxe is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A (not applicable for 120v rated devices)	
Voltage fluctuations/flicker emissions IEC61000-3-3	Complies (applicable to 50Hz rated devices only)	

4.

Guidance and manufacturer's declaration – electromagnetic immunity				
The Innospire Deluxe is intended for use in the electromagnetic environment specified below. The customer or the user of the Innospire Deluxe should assure that it is used in such an environment.				
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance	
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.	
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines and dc lines	±2 kV for power supply lines and dc lines	Mains power quality should be that of a typical commercial or hospital environment.	
Surge IEC 61000-4-5	±1 kV line(s) and neutral ±2 kV line(s) to earth	±1 kV line(s) and neutral ±2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.	
Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11	<5 % U _T (>95 % dip in U _T) for 0.5 cycle 40 % U _T (60 % dip in U _T) for 5 cycles 70 % U _T (30% dip in U _T) for 25 cycles <5 % U _T (>95 % dip in U _T) for 5s	<5 % U _T (>95 % dip in U _T) for 0.5 cycle 40 % U _T (60 % dip in U _T) for 5 cycles 70 % U _T (30% dip in U _T) for 25 cycles <5 % U _T (>95 % dip in U _T) for 5s	Mains power quality should be that of a typical commercial or hospital environment. If a dips or an interruption of mains power occurs, the current of the Innospire Deluxe may be dropped off from normal level, it may be necessary to use uninterruptible power supply or a battery.	
Power frequency (50/60-Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.	

NOTE U_T is the a.c. mains voltage prior to application of the test level.

5.

Guidance and manufacturer's declaration – electromagnetic immunity			
The Innospire Deluxe is intended for use in the electromagnetic environment specified below. The customer or the user of the Innospire Deluxe should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3Vrms 150kHz to 80MHz	3Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the Innospire Deluxe, including cables, than the recommended separation distance, calculated from the equation applicable to the frequency of the transmitter: Recommended separation distance d = 1.2√P d=1.2√P 80MHz to 800MHz d=2.3√P 800MHz to 2.5GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3V/m 80MHz to 2.5GHz	3V/m	
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Innospire Deluxe is used exceeds the applicable RF compliance level above, the Innospire Deluxe should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Innospire Deluxe. b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

192

6.

Recommended separation distances between portable and mobile RF communications equipment and the Innospire Deluxe			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz d=1.2√P	80 MHz to 800 MHz d=1.2√P	800 MHz to 2.5 GHz d=2.3√P
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

193

InnoSpire Deluxe is a trademark of Koninklijke Philips Electronics N.V. and its affiliates, all rights reserved.

© 2013 Koninklijke Philips Electronics N.V. All rights are reserved.

CE 0123

Manufactured for

Respironics Respiratory Drug Delivery (UK) Ltd
Chichester Business Park, City Fields Way
Tangmere, Chichester, PO20 2FT UK

Made in China

www.philips.com/respironics

PN 1113054 R00