# Snabbguide MEMOplanner

---

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Sida 3 (8)

5. Välj datum.
6. Välj typ.

The image illustrates the steps for adding an activity in a mobile application. It consists of three sequential screenshots:

- Calendar Selection:** A calendar for October 2016 is shown. The dates from the 1st to the 30th are highlighted in pink, indicating they are available for selection. The interface includes a back arrow, a title '+ Ny aktivitet Välj datum', and a forward arrow.
- Activity Form:** A form titled '+ Ny aktivitet' is displayed. It has tabs for 'Namn', 'Extra', 'Larm', and 'Påminn'. The form contains several fields: 'Namn', 'Bild', 'Datum' (set to 'tisdag 27 september'), 'Typ' (set to 'Höger'), 'Tid', and two checkboxes: 'Kvitterbar' and 'Ta bort efteråt'. Navigation arrows are present for the first four fields. At the bottom, there are '< Föregående' and 'Nästa >' buttons.
- Type Selection:** A screen titled '+ Välj typ' shows three radio button options: 'Heldag', 'Vänster', and 'Höger'. The 'Höger' option is selected. At the bottom, there are 'Avbryt' and 'OK' buttons.

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

The image shows a mobile application interface for adding a new activity. The main screen is titled 'Ny aktivitet Återkommande' and has four radio button options: 'Enstaka', 'Veckovis', 'Varannan vecka', and 'Månadsvis'. The 'Veckovis' option is selected. Below these options are 'Månadsvis' and 'Årsvis' options. At the bottom, there is a green 'Spara' button circled in red. To the left, a 'Månadsvis' calendar view is shown with a grid of days from 1 to 31. To the right, a 'Veckovis' weekly view is shown with days of the week (mån, tis, ons, tor, fre, lör, sön) and a 'Markera alla' button. Arrows point from the 'Veckovis' option on the main screen to the weekly view, and from the 'Spara' button to the main screen.

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