Guardian™ Connect

Transmitter

Sender

Lähetin

Sändare

Sender

Nadajnik

Contacts:

Africa: Medtronic Africa (Pty) Ltd.

Tel: +27 (0) 11 677 4800

Argentina: Corpomedica S.A.

Tel: +(11) 4 814 1333

Medtronic Directo 24/7: +0800 333 0752

Australia: Medtronic Australasia Pty. Ltd.

Tel: 1800 668 670 (product orders)

Tel: 1800 777 808 (customer help)

Azerbaijan: Isomed

Tel: +994 (12) 464 11 30

Bangladesh: Sonargaon Healthcare Pvt Ltd.

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Fax: +34 91 625 03 90

24 horas: +34 900 120 330

Europe: Medtronic Europe S.A. Europe, Middle
East and Africa Headquarters

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Malaysia: Medtronic International Ltd.
Tel: +603 7946 9000

Middle East and North Africa: Regional Office
Tel: +961-1-370 670

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Ticaret Ltd. Sirketi.

Tel: +90 216 4694330

USA: Medtronic Diabetes Global Headquarters

24 Hour HelpLine: +1-800-646-4633

To order supplies: +1-800-843-6687

Ukraine: Medtronic B.V. Representative office in
Ukraine

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Non-stop helpLine:

+420 233 059 950

The Guardian™ Connect transmitter is a component of the Guardian™ Connect continuous glucose monitoring (CGM) system. The transmitter is compatible only with the Enlite™ glucose sensor. The transmitter receives data from the sensor. The transmitter also receives and sends data to the Guardian™ Connect app through a Bluetooth® Smart wireless connection.

transmitter

tester

charger

Guardian Connect transmitter kit components

A complete transmitter kit includes the following components:

- Guardian Connect transmitter (MMT-7821)
- Two testers (MMT-7736)
- Charger (MMT-7715)

Indications for use

The transmitter is indicated for single-patient or multiple-patient use as a component of select Medtronic CGM systems.

Contraindications

Do not expose your transmitter to MRI equipment, diathermy devices, or other devices that generate strong magnetic fields. If your transmitter is inadvertently exposed to a strong magnetic field, discontinue use and contact the 24 Hour HelpLine or your local representative for further assistance.

Warnings

- Always refer to the sensor user guide for all precautions, warnings, and instructions relating to the sensor. Not referring to the sensor user guide can result in serious injury or damage to the sensor.
- Do not allow children to put small parts in their mouth. This product poses a choking hazard for young children.

- Do not change or modify the device unless expressly approved by Medtronic Diabetes. Modifying the device can cause serious injury, interfere with your ability to operate the device, and void your warranty.
- Do not use the tester if it comes in contact with blood. Touching blood can cause infection. Dispose of the tester according to the local regulations for medical waste disposal, or contact your healthcare professional for disposal information.
- Bleeding may occur after inserting the sensor. Always make sure that the site is not bleeding before connecting the transmitter to the sensor. Blood can get into the transmitter connector and damage the device. Discard the device if damaged. If bleeding occurs, apply steady pressure with a sterile gauze or clean cloth at the insertion site until bleeding stops. After bleeding stops, connect the transmitter to the sensor.
- Contact the 24 Hour HelpLine or your local representative if you experience any adverse reactions associated with the transmitter or sensor. Adverse reactions can cause serious injury.

Precautions

- Only use the Enlite (MMT-7008) glucose sensor with the transmitter. Do not use any other sensor. Other sensors are not intended for use with the transmitter.
- Only use the tester (MMT-7736) with the transmitter. Do not use any other test plug. Other test plugs are not intended for use with the transmitter.
- Always use the tester when cleaning the transmitter. Do not use any other test plug with the transmitter. Use of another test plug can allow water to get into the transmitter. Water can damage the transmitter.
- Do not twist the tester or sensor while attached to the transmitter. Twisting the tester or sensor will damage the transmitter.
- Do not allow the tester to come in contact with any liquid when not connected to the transmitter. A wet tester can damage the transmitter.
- Do not allow the transmitter to come in contact with any liquid when not connected to a sensor or to the tester. Moisture will damage the transmitter and a wet transmitter can damage the sensor.
- Do not clean the o-rings on the tester with any substances. Cleaning the o-rings can damage the tester.

Exposure to magnetic fields and radiation

- Do not expose your transmitter to Magnetic Resonance Imaging (MRI) equipment, diathermy devices, or other devices that generate strong magnetic fields (for example, x-ray, CT scan, or other types of radiation). The strong magnetic fields can

cause the device to malfunction, and result in serious injury. If your transmitter is exposed to a strong magnetic field, discontinue use and contact the 24 Hour HelpLine or your local representative for further assistance.

- Always remove your sensor and transmitter before entering a room that has x-ray, MRI, diathermy, or CT scan equipment. This equipment has strong magnetic fields that can cause the device to malfunction, and result in serious injury. If your sensor or transmitter is exposed to a strong magnetic field, discontinue use and contact the 24 Hour HelpLine or your local representative for further assistance.

Radio Frequency (RF) communication

This device complies with Part 15 of the FCC Rules. Operation is subject to two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

This device has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This device generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this device does cause harmful interference to radio or television reception, which can be determined by turning the device off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the device and the receiver.
- Decrease the distance between the transmitter and the compatible mobile device to 6.1 meters (20 feet) or less.
- Increase the separation between the transmitter and the equipment that is receiving/emitting interference.

Do not change or modify the internal RF transmitter or antenna unless expressly approved by Medtronic Diabetes. Doing so could interfere with your ability to operate the device.

Directive 1999/5/EC

Medtronic declares that this product is in conformity with the essential requirements of Directive 1999/5/EC on Radio and Telecommunications Terminal Equipment.

For additional information, contact Medtronic MiniMed at the address or phone number provided on the back cover.

IEC60601-1-2; Special EMC Precautions for Medical Electrical Equipment

- 1 Special Precautions regarding Electromagnetic Compatibility (EMC): This body worn device is intended to be operated within a reasonable residential, domestic, public or work environment, where common levels of radiated “E” (V/m) or “H” fields (A/m) exist; such as cellular phones, WiFi, Bluetooth, electric can openers, microwave and induction ovens. This device generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the provided instructions, may cause harmful interference to radio communications.
- 2 Portable and mobile RF communications equipment can affect Medical Electrical Equipment communications as well. If you encounter RF interference from a mobile or stationary RF transmitter, move away from the RF transmitter that is emitting the interference.

Assistance

Please contact your local representative for assistance. Refer to the Medtronic Diabetes International Contacts list at the beginning of this user guide for contact information.

Charger

The transmitter contains a non-replaceable, rechargeable battery that you can recharge as needed with the charger. The charger has a green light that shows the charging status and a red light that communicates any problems during charging. If you see a red light, see the Troubleshooting section. The charger needs one AAA alkaline battery.


Note: *If the battery is installed incorrectly or is low, the charger will not work. Repeat the battery installation steps using a new battery.*

Installing a battery in the charger

To install a battery in the charger:

- 1 Push the battery cover in and slide it off (as shown in the image in step 3).
- 2 Insert a new AAA alkaline battery. Make sure the + and - symbols on the battery align with these same symbols shown on the charger.

- Slide the cover back on the charger until it clicks into place.



Charging the transmitter

Caution: Always charge the transmitter before inserting your sensor. A depleted transmitter does not function. A fully charged transmitter works at least six days without recharging. A depleted transmitter can take up to two hours to recharge.

Note: Do not store the transmitter on the charger for more than 60 days. Disconnect and reconnect to the charger to re-charge again before use. If the transmitter is left on the charger for more than 60 days, the battery will be permanently damaged.

To charge the transmitter:

- Push the two components together to connect the transmitter to the charger.
- Within 10 seconds after the transmitter is connected, a green light on the charger will flash for one to two seconds as the charger powers on. For the rest of the charging time, the green light on the charger will continue to flash in a pattern of four flashes with a pause between the four flashes.
- When charging is complete, the green light on the charger will stay on, without flashing, for 15 to 20 seconds and then turn off.
- After the green charger light turns off, disconnect the transmitter from the charger. The green light on the transmitter will flash 10 times and then turn off.

Pairing your transmitter

Always refer to the system user guide for instructions on how to pair your transmitter to your mobile device.

Inserting the sensor


Always refer to yourserter user guide for instructions on how to insert the sensor.

Connecting the transmitter to the sensor

Before proceeding, have your system user guide available.

To connect the transmitter to the sensor:

- 1 After the sensor is inserted, consult yourserter user guide for details on applying the required tape before connecting the transmitter.
- 2 Hold the rounded end of the inserted sensor to prevent it from moving during connection.
- 3 Hold the transmitter as shown. Line up the two notches on the transmitter with the side arms of the sensor. The flat side of the transmitter should face the skin.
- 4 Slide the transmitter onto the sensor connector until the sensor arms snap into the notches on the transmitter. If the transmitter is properly connected, and if the sensor has had enough time to become hydrated, the green light on the transmitter will flash 6 times.

Note: If the transmitter does not flash, see *Troubleshooting*, on page 19.


- 5 When the transmitter light flashes green after connecting to the sensor, follow the prompts on the app to select new or existing sensor. For more instructions, see your system user guide.
- 6 Attach the adhesive tab of the sensor to the transmitter.
- 7 Refer to theserter user guide for instructions on how to apply additional tape if necessary.
- 8 Follow the instructions that appear on the app or follow the instructions in your system user guide.

Disconnecting the transmitter from the sensor

Before proceeding, have your system user guide available.

To disconnect the transmitter from the sensor:

- 1 Carefully remove any tape from the transmitter and sensor.
- 2 Remove the adhesive tab from the top of the transmitter.
- 3 Hold the transmitter as shown, and pinch the flexible side arms of the sensor between your thumb and forefinger.
- 4 Gently pull the transmitter away from the sensor.
- 5 Follow the instructions that appear in the app or follow the instructions in your system user guide.

Bathing and swimming

After the transmitter and sensor are connected, they form a waterproof seal to a depth of 2.4 meters (eight feet) for up to 30 minutes. You can shower and swim without removing them.


Tester

The tester is used to test the transmitter to make sure it is working. It is also used as a required component for cleaning the transmitter. Properly connecting the tester to the transmitter ensures that fluids do not come in contact with the connector pins inside the transmitter. Fluids can cause connector pins to corrode and affect the performance of the transmitter.


Do not twist the tester while attached to the transmitter. This will damage the transmitter.

The tester can be used for one year. If you continue to use the tester for more than one year, the connector pins inside the transmitter could be damaged, because the tester cannot continue to provide a waterproof seal. For instructions on how to check the connector pins, see *Inspecting the transmitter connector pins*, on page 8.

Caution: Only use the tester (MMT-7736) with the transmitter. Do not use any other test plug. Other test plugs are not intended for use with the transmitter, and will damage the transmitter and the tester.


Correct tester has gold on both sides of connector.

Inspecting the transmitter connector pins

This image is an example of how the connector pins should look.

Look inside the connector opening of the transmitter to make sure that the connector pins are not damaged or corroded. If the connector pins are damaged or corroded, the transmitter cannot communicate with the charger, monitor, or pump. Contact the 24 Hour HelpLine or your local representative. It may be time to replace your transmitter.

Also look for moisture inside the connector opening. If you see any moisture, allow the transmitter to dry for at least one hour. Moisture inside the connector opening could cause the transmitter to not work properly, and could cause corrosion and damage over time.


Connecting the tester for testing or cleaning

Before proceeding, have your system user guide available.

To connect the tester:

- 1 Hold the transmitter and the tester as shown. Line up the flat side of the tester with the flat side of the transmitter.
- 2 Push the tester into the transmitter until the flexible side arms of the tester click into the notches on both sides of the transmitter.

When properly connected, the green light on the transmitter flashes 6 times.


- 3 To test the transmitter, check the sensor icon in the app to ensure that the transmitter is sending a signal (see your system user guide).
- 4 To clean the transmitter, see *Cleaning the transmitter*, on page 9.
- 5 After testing or cleaning, disconnect the tester from the transmitter.

Disconnecting the tester

To disconnect the tester:

- 1 Hold the transmitter body as shown and pinch the side arms of the tester.
- 2 With the tester arms pinched, gently pull the transmitter away from the tester.

Note: To save transmitter battery life, do NOT leave the tester connected after cleaning or testing.

Cleaning the transmitter

The transmitter is intended for personal use at home (single-patient use) or for use in healthcare facilities (multiple-patient use). Single-patient use requires cleaning after each use, while multiple-patient use requires cleaning and disinfection after each use. When using the transmitter in a healthcare facility, always follow the cleaning and disinfecting procedure for multiple-patient use.

Warning: Do not discard the transmitter in a medical waste container or expose it to extreme heat. The transmitter contains a battery that may ignite, and result in serious injury.

Note: The tester is a required component for cleaning the transmitter. For details, see *Tester*, on page 7.

For single-patient use

Always clean the transmitter after each use.

To clean the transmitter, you need the following materials:

- mild liquid soap
- soft-bristled toddler toothbrush
- container
- clean and lint-free dry cloths


Use life

The transmitter can be cleaned up to 122 cycles or one year, whichever comes first. Discard the transmitter at this point. If you continue to use the transmitter beyond 122 cycles or one year, the cleaning process may damage the device. Contact Medtronic to order a new transmitter.

Warning: Do not use the device if you see any cracking, flaking, or damage to the housing. Cracking, flaking, or damage to the housing are signs of deterioration. Deterioration of the housing can affect the ability to properly clean the transmitter, and result in serious injury. Call the 24 Hour HelpLine or your local representative and discard the device according to local regulations for battery disposal (nonincineration), or contact your healthcare professional for disposal information.


To clean the transmitter:

- 1 Wash your hands thoroughly.
- 2 Attach the tester to the transmitter.


- 3 If there is adhesive residue on the transmitter, see *Removing adhesive residue*, on page 17.

- 4 Rinse the transmitter under room temperature tap water for at least one minute, and until visibly clean. Make sure all hard-to-reach areas are rinsed completely.


- 5 Prepare a mild liquid soap solution using 5 milliliters (1 teaspoon) of mild liquid soap per 3.8 liters (1 gallon) of room temperature tap water.
- 6 With the tester still attached, submerge the transmitter in the cleaning solution and soak for one minute.


- 7 Holding the tester, brush the entire surface of the transmitter using a soft-bristled toddler toothbrush. Make sure to brush all hard-to-reach areas until visibly clean.

- 8 Rinse the transmitter under running room temperature tap water for at least one minute, and until all visible liquid soap is gone.


- 9 Dry the transmitter and tester with a clean, dry cloth.

- 10 Place the transmitter and tester on a clean, dry cloth and air dry them completely.

- 11 Disconnect the tester from the transmitter by gently squeezing the arms of the tester.

For multiple-patient use

When using the transmitter in a healthcare facility, always clean and disinfect the transmitter after each use.

Warning: You must adhere to Standard Precautions when handling or using this device. All parts of the system should be considered potentially infectious and are capable of transmitting blood-borne pathogens between patients and healthcare professionals.

The transmitter must be disinfected after use on each patient. This system may only be used for testing multiple patients when Standard Precautions and Medtronic's disinfection procedures are followed.

To clean the transmitter, you need the following materials:

- gloves
- mild liquid soap
- soft-bristled toddler toothbrush
- 8% bleach
- two containers
- clean and lint-free dry cloths

Use life


The transmitter can be cleaned and disinfected up to 122 cycles or one year, whichever comes first. Discard the transmitter at this point. If you continue to use the transmitter beyond 122 cycles or one year, the cleaning and disinfection process may damage the device. Contact Medtronic to order a new transmitter.

To clean and disinfect the transmitter:


- 1 Wash your hands and put on gloves.
- 2 Inspect the inside of the transmitter's connector opening for any sign of body fluid. For instructions on how to inspect the connector pins, see *Inspecting the transmitter connector pins*, on page 8.

Caution: The person inspecting the transmitter must have sufficient vision that enables him or her to see small drops of body fluid or debris.


Warning: If you see any body fluid in the connector opening, you must discard the transmitter. Because the transmitter contains a battery, do not discard in a bio-waste container. Instead, continue to clean and disinfect the transmitter, and then discard according to local regulations for battery disposal (non-incineration).


3 Attach the tester to the transmitter.

- 4 If there is adhesive residue on the transmitter, see *Removing adhesive residue*, on page 17.
- 5 Rinse the transmitter under room temperature tap water for at least one minute and until visibly clean. Make sure all hard-to-reach areas are rinsed completely.


- 6 Prepare a mild liquid soap solution using 5 milliliters (1 teaspoon) of mild liquid soap per 3.8 liters (1 gallon) of room temperature tap water. Make sure to prepare a fresh solution for each use.
- 7 With the tester still attached, fully submerge the transmitter in the cleaning solution and soak for one minute.


- 8 Holding the tester, brush the entire surface of the transmitter using a soft-bristled toddler toothbrush. Make sure to brush all hard-to-reach areas until visibly clean.


- 9 Rinse the transmitter under running room temperature tap water for at least one minute, and until all visible liquid soap is gone.


- 10 Dry the transmitter and tester with a clean, dry cloth.

- 11 Prepare a 1:10 bleach solution by using one (1) part 8.25% bleach to nine (9) parts water, for a final concentration of 0.8%. Make sure to prepare a fresh solution for each use.
- 12 Make sure that you have completed the previous cleaning steps before disinfection. With the tester still attached, soak the transmitter in the bleach solution for 20 minutes.


- 13 Rinse the transmitter under running room temperature tap water for three minutes.

- 14 Place the transmitter and tester on a clean, dry cloth and air dry them completely.

Warning: If you saw any body fluid inside the connector opening on earlier inspection, you must now discard the transmitter with tester still attached, according to local regulations for battery disposal (non-incineration).

- 15 Disconnect the tester from the transmitter by gently squeezing the arms of the tester.

- 16 Inspect the housing of the transmitter for any signs of cracking, flaking, or damage. If you see any of these signs, you must now discard the disinfected transmitter according to local regulations for battery disposal (non-incineration).

Warning: Do not use the device if you see any cracking, flaking, or damage to the housing. Cracking, flaking, or damage to the housing are signs of deterioration. Deterioration of the housing can affect the ability to properly clean the transmitter, and result in serious injury. Call the 24 Hour HelpLine or your local representative and discard the device according to local regulations for battery disposal (nonincineration), or contact your healthcare professional for disposal information.

- 17 Discard the used gloves and thoroughly wash hands with soap and water.

Removing adhesive residue


You may need to perform this procedure if there is adhesive residue present on the transmitter. If you visually inspect the transmitter and see adhesive residue on it, follow the instructions below.

To remove adhesive residue, you need medical adhesive remover (such as Detachol® - which is a mineral spirit) and cotton swabs.

Note: During testing, Medtronic MiniMed used Detachol to remove the adhesive residue from the transmitter. Detachol is recommended for use but may not be available in all countries.

To remove adhesive residue:

- 1 Make sure the tester is attached to the transmitter.
- 2 Holding the tester, saturate a cotton swab in the Detachol solution and gently rub the adhesive residue on the transmitter until it is fully removed.

- 3 Continue with the appropriate single-patient or multiple-patient cleaning procedure. See *Cleaning the transmitter*, on page 9 for details.

Cleaning the charger

This procedure is for general cleaning as required, based on physical appearance.

Caution: Do not immerse the charger in water or any other cleaning agent. The charger is not waterproof. Water can damage the charger, and cause the device to malfunction.

Warning: Dispose the charger according to the local regulations for battery disposal, or contact your healthcare professional for disposal information. The charger may ignite upon incineration.

Warning: For multiple-patient use, always clean and disinfect the transmitter after removing it from the patient and before attaching it to the charger. If blood comes in contact with any surface areas of the charger, the contaminated device must be discarded. The charger contains a battery which may explode upon incineration.

To clean the charger:

- 1 Wash your hands thoroughly.

- 2 Use a damp cloth with mild cleaning solution, such as a dishwashing detergent, to clean any dirt or foreign material from the outside of the charger. Never use organic solvents, such as paint thinner or acetone, to clean the charger.
- 3 Place the charger on a clean, dry cloth and air dry for two to three minutes.

Troubleshooting

The following table contains troubleshooting information for the transmitter, charger, and tester. For more information about troubleshooting, see your system user guide.

Problem	Likely Cause(s)	Resolution
You connected the transmitter to the charger and no lights came on.	The transmitter connector pins are damaged or corroded. Your charger battery has no power or no battery is inserted.	1 Check the transmitter connector pins for damage or corrosion. For more information about your connector pins, see <i>Inspecting the transmitter connector pins</i>, on page 8. If the pins are damaged or corroded, contact the 24 Hour HelpLine or your local representative. It may be time to replace your transmitter. 2 If there is no damage to the connector pins, replace the battery in the charger. For instructions on replacing your charger battery, see <i>Installing a battery in the charger</i>, on page 4.
During charging, the flashing green light on the charger turns off and you see a longer flashing red light on the charger.	Your charger battery is low on power.	Replace the battery in the charger. For instructions on replacing your charger battery, see <i>Installing a battery in the charger</i> , on page 4.
During charging, the flashing green light on the charger turns off and you see a series of quick flashing red lights on the charger for two seconds at a time.	Your transmitter is low on power.	1 Charge the transmitter continuously for one hour. If flashing does not stop, proceed to step 2. 2 Charge the transmitter continuously for eight hours. If flashing does not stop, call the 24 Hour HelpLine or your local representative. It may be time to replace your transmitter.

Problem	Likely Cause(s)	Resolution
During charging, a mix of quick and long flashing red lights appear on the charger.	Your charger <i>and</i> your transmitter are low on power.	1 Replace the battery in the charger. For instructions on replacing your charger battery, see <i>Installing a battery in the charger, on page 4</i>. 2 Charge the transmitter continuously for one hour. If the quick flashing red lights do not stop, proceed to step 3. 3 Charge the transmitter continuously for eight hours. If flashing does not stop, call the 24 Hour HelpLine or your local representative. It may be time to replace your transmitter.
The green light on the transmitter does not flash when you connect it to the sensor.	Your transmitter is not fully connected. Your transmitter is low on power. Your sensor is not properly inserted into your body.	1 Disconnect the transmitter from the sensor. 2 Wait for five seconds and reconnect them. If the green light still does not flash, proceed to step 3. 3 Fully charge the transmitter and connect it to the tester. If the green light still does not flash, see troubleshooting on “The green light on the transmitter does not flash when you connect it to the tester”. If the green light flashes, proceed to step 4. 4 Disconnect the transmitter from the tester, wait at least five seconds, and connect the transmitter to the sensor. If the green light still does not flash, proceed to step 5. 5 The sensor may not be properly inserted into your body. Remove the sensor from your body and insert a new sensor.
The green light on the transmitter does not flash when you connect it to the tester.	Your transmitter is not fully connected. Your transmitter is low on power.	1 Check the connection between the transmitter and the tester. If the green light still does not flash, proceed to step 2. 2 Fully charge the transmitter. 3 Test the transmitter with the tester again. If you still do not see the green light flash, call the 24 Hour HelpLine or your local representative. It may be time to replace your transmitter.

Problem	Likely Cause(s)	Resolution
Your transmitter battery does not last for six days.	Your transmitter is not fully charged when you connect it to the sensor. The transmitter and the app on your compatible mobile device frequently lose wireless connection.	1 Fully charge the transmitter before connecting it to the sensor. If the transmitter battery still does not last for six days, proceed to step 2. 2 Move away from any equipment that can cause RF interference. For more information on RF interference, see <i>Radio Frequency (RF) communication, on page 3</i>. 3 Make sure your compatible mobile device and your transmitter are located on the same side of your body to minimize any RF interference. If your fully charged transmitter battery continues to lose power before a full six days, call the 24 Hour HelpLine or your local representative. It may be time to replace your transmitter.
Your transmitter has lost connection with the app on your compatible mobile device. Note: <i>An alert occurs and a message appears when your transmitter has lost connection with the app on your compatible mobile device.</i>	Your compatible mobile device is out of range. There is RF interference from other devices.	1 Move away from any equipment that can cause RF interference. For more information on RF interference, see <i>Radio Frequency (RF) communication, on page 3</i>. If your transmitter is still not communicating with the app on your compatible mobile device, proceed to step 2. 2 Make sure your compatible mobile device and your transmitter are located on the same side of your body to minimize any RF interference. If your transmitter is still not communicating with the app, call the 24 Hour HelpLine or your local representative for assistance.

Storing the devices

Store the transmitter, charger, and tester in a clean, dry location at room temperature. If the transmitter is not in use, you must charge the transmitter at least once every 60 days. Although not required, you may store the transmitter on the charger. If you are storing the transmitter on the charger, you must disconnect and reconnect the charger and the transmitter at least once every 60 days. If the transmitter is left on the charger for more than 60 days, the battery will be permanently damaged.

Disposal

Discard the transmitter according to local regulations for battery disposal, or contact your healthcare professional for disposal information.

Specifications

Biocompatibility	Transmitter: Complies with EN ISO 10993-1
Applied parts	Transmitter Sensor
Operating conditions	Transmitter temperature: 0 °C to 45 °C (32 °F to 113 °F) Caution: When operating the transmitter on a tester in air temperatures greater than 41 °C (106 °F), the temperature of the transmitter may exceed 43 °C (109 °F) Transmitter relative humidity: 10% to 95% with no condensation Transmitter pressure: 57.60 to 106.17 kPa (8.4 to 15.4 psi) Charger temperature: 10 °C to 40 °C (50 °F to 104 °F) Charger relative humidity: 30% to 75% with no condensation
Storage conditions	Transmitter temperature: -20 °C to 55 °C (-4 °F to 131 °F) Transmitter relative humidity: up to 95% with no condensation Transmitter pressure: 57.60 to 106.17 kPa (8.4 to 15.4 psi) Charger temperature: -10 °C to 50 °C (14 °F to 122 °F) Charger relative humidity: 10% to 95% with no condensation
Battery life	Transmitter: Six days of continuous glucose monitoring immediately following a full charge Charger: The charger uses one new AAA battery to charge the transmitter
Transmitter frequency	2.4 GHz band, Bluetooth Smart (version 4.0)
Maximum output power (EIRP)	0.1 mW (-9.9 dBm)
Operating range	Up to 6.1 meters (20 feet)
Transmitter expected service life	The transmitter expected service life is one year depending on patient usage.

Transmitter wireless communication

Quality of service

The transmitter and app connect via a BLE network. The transmitter sends glucose data and system related alerts to the app. The transmitter and the app verify the integrity of received data after wireless transmission. Quality of the connection is in accordance with the Bluetooth Specification v4.0.

Data security

The transmitter is designed to only accept radio frequency (RF) communications from recognized and linked devices. You must program the app to accept information from a specific transmitter. Transmitted sensitive data is encrypted to prevent unauthorized receipt or communication.

Guidance and manufacturer's declaration

Guidance and Manufacturer's Declaration - Electromagnetic Emissions		
The transmitter is intended for use in the electromagnetic environment specified below. The customer or the user of the transmitter should make sure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic Environment - Guidance
RF emissions CISPR 11	Group 1	The transmitter must emit electromagnetic energy in order to perform its intended function. Nearby electronic equipment may be affected.
RF emissions CISPR 11	Class B	The transmitter is suitable for use in all establishments, including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

Guidance and Manufacturer's Declaration - Electromagnetic Immunity			
The transmitter is intended for use in the electromagnetic environment specified below. The customer or the user of the transmitter should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±2 kV, ±4 kV, ±8 kV Air ±2 kV, ±4 kV, ±6 kV Indirect	±2 kV, ±4 kV, ±8 kV Air ±2 kV, ±4 kV, ±6 kV Indirect	For use in a typical domestic, commercial, or hospital environment.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	Not applicable	Requirement does not apply to this battery powered device.
Surge IEC 61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	Not applicable	Requirement does not apply to this battery powered device.

Guidance and Manufacturer's Declaration - Electromagnetic Immunity			
The transmitter is intended for use in the electromagnetic environment specified below. The customer or the user of the transmitter should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Voltage dips, short interruptions and voltage variations on power supply lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycle	Not applicable	Requirement does not apply to this battery powered device.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	400 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical domestic, commercial, or hospital environment.
Note: U_T is the a.c. mains voltage prior to application of the test level.			

Guidance and Manufacturer's Declaration - Electromagnetic Immunity			
The transmitter is intended for use in the electromagnetic environment specified below. The customer or user of the transmitter should assure that it is used in such an electromagnetic environment.			
Immunity Test	IEC 60601 Level	Compliance Level	Electromagnetic Environment Guidance
Conducted RF IEC 61000-4-6	3 V/m 150 kHz to 80 MHz	Not applicable	Not applicable

Guidance and Manufacturer's Declaration - Electromagnetic Immunity

The transmitter is intended for use in the electromagnetic environment specified below. The customer or user of the transmitter should assure that it is used in such an electromagnetic environment.

Immunity Test	IEC 60601 Level	Compliance Level	Electromagnetic Environment Guidance
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 6 GHz	10 V/m 80 MHz to 6 GHz	Portable and mobile RF communications equipment should be used no closer to any part of the transmitter, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Refer to the recommended separation distance table for more information. $d = 0.35\sqrt{P}$ 80 MHz to 800 MHz $d = 0.70\sqrt{P}$ 800 MHz to 6 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b. Interference may occur in the vicinity of equipment marked with the following symbol: 

Note: At 80 MHz and 800 MHz, the higher frequency range applies.

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption, and reflection from structures, objects and people.

Guidance and Manufacturer's Declaration - Electromagnetic Immunity

The transmitter is intended for use in the electromagnetic environment specified below. The customer or user of the transmitter should assure that it is used in such an electromagnetic environment.

Immunity Test	IEC 60601 Level	Compliance Level	Electromagnetic Environment Guidance
---------------	-----------------	------------------	--------------------------------------

^aField strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcasts and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the transmitter is used exceeds the applicable RF compliance level above, the transmitter should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the transmitter.

^bOver the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the transmitter

The transmitter is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the transmitter users can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment and the transmitter as recommended below, according to the maximum output power of the communications equipment.


Rated maximum output power of transmitter (W)	Separation distance according to the frequency of transmitter (m)		
	150 kHz to 80 MHz Not applicable	80MHz to 800MHz $d = 0.35 \sqrt{P}$	800MHz to 6.0GHz $d = 0.70 \sqrt{P}$
0.01	Not applicable	0.035	0.07
0.1	Not applicable	0.11	0.22
1	Not applicable	0.35	0.7
10	Not applicable	1.1	2.2
100	Not applicable	3.5	7

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

Icon Table

	Serial number
	Catalogue or model number
(1x)	One transmitter, charger, and senter per container/package
(2x)	Two testers per container/package
	Date of manufacture (YYYY-MM-DD)
	Manufacturer
	Must refer to instruction manual before every use (appears blue on label).
	Temperature limit
	Non-ionizing electromagnetic radiation (RF communication)
	Configuration or unique version identifier
	Degree of protection against electric shock: Type BF applied part
IP48	Transmitter: 4 is the level of protection against solid objects with a diameter above 1mm. 8 is the level of protection against the effects of continuous immersion in water [2.4 meters (8 feet) immersion for 30 minutes].
	Humidity limitation
	This product conforms to Australia Radio Requirements
CE 0459	Marking of Conformity: This symbol means the device fully complies with MDD 93/42/EEC (NB 0459).
	Authorized representative in the European community
	Fragile, handle with care
	Keep dry

	Recycle cardboard, paper, plastic packaging supplies and unwanted written material.
 	WEEE Initiative: DO NOT THROW IN TRASH. Recycle device according to local disposal requirements.
	Magnetic Resonance (MR) unsafe: keep away from magnetic resonance imaging (MRI) equipment
	<i>Bluetooth</i> ® wireless technology or <i>Bluetooth</i> ® enabled


©2016 Medtronic MiniMed, Inc. All rights reserved.

MiniMed™, Enlite™, and Guardian™ are trademarks of Medtronic MiniMed, Inc.

Bluetooth® is a registered trademark of Bluetooth SIG, Inc.

Detachol® is a registered trademark of Ferndale Laboratories Inc.

Guardian™ Connect-senderen er en del av Guardian™ Connect-systemet for kontinuerlig glukosemåling (CGM). Senderen er kun kompatibel med Enlite™ - glukosesensoren. Senderen mottar data fra sensoren. Senderen mottar og sender også data til Guardian™ Connect-appen via en Bluetooth® Smart trådløs tilkobling.

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©2016 Medtronic MiniMed, Inc. Wszelkie prawa zastrzeżone.
MiniMed™, Enlite™ i Guardian™ są znakami towarowymi firmy Medtronic MiniMed, Inc.
Bluetooth® jest zastrzeżonym znakiem towarowym firmy Bluetooth SIG, Inc.
Detachol® jest zastrzeżonym znakiem towarowym firmy Ferndale Laboratories Inc.