TS-150 Time Delay

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
Contents

User manual for the TS-150	30
What's in the box	30
Daily use	30
Use of microphone	30
Lamps and controls	31
Interference sources and troubleshooting	31
Troubleshooting table	32
AGC (Automatic gain control)	32
Supplier	33
Installation instructions for the TS-150	33
Preparation	33
Mounting	33
Connections	34
Adjustments	35
Adjusting the input levels	35
Microphone for the TS-150	37
Adjusting the sound	37
AGC or not?	38
Adjusting the loop current (GAIN)	38
Selection of loop cable	39
Notes	39
Technical data	40
Cleaning, safety and maintenance	41
Drilling template	42

User manual for the TS-150

What's in the box

Daily use

The TS-150 has no ON/OFF switch and may be continuously connected to the power supply. After installation, it is no longer necessary to make any adjustments to the amplifier. Sound sources such as TV sets, radios and microphones are merely switched on and off. The hearing instrument is set to the T or MT position.

Only one sound source plus the microphone should be on at a time.


Synchronization of sound - the DELAY knob

Unsynchronized sound can be found in certain LCD or PLASMA TV sets. If the sound is not synchronized with the image (the sound precedes the image) you can adjust this with the knob marked "DELAY". Increase the delay of the sound by turning the DELAY knob clockwise until you find that the sound is synchronized with the image. This should only be needed to do once.

Use of microphone

A microphone can be connected to the TS-150 via the microphone socket. Hold the microphone close to your mouth when speaking. When you talk into the microphone, the signals from the other sound sources are greatly weakened so the sound from the microphone dominates. The sound from the other sources returns to its full volume a few seconds after you have stopped talking into the microphone.

Lamps and controls

Interference sources and troubleshooting

Many interference sources can affect the sound quality of an inductive loop system. As the sound is transmitted wirelessly via electromagnetic waves, electrical equipment in the vicinity of the receiver (hearing instrument) can cause buzzing and interference. Typical interference sources include heat cables, fluorescent tubes, computer screens and TV sets. Before installing an inductive loop, you should first check for any interference sources.

If a problem occurs with the normal use of the system, you should check that the lamps on the amplifier light up as described above, that the receiver or hearing instrument is switched on and that its battery is operational. Consult the table on the next page before calling out your supplier.


If you need further help with troubleshooting, you should contact your supplier!

Troubleshooting table

Fault	Possible cause	Action
Sound is faint.	1. Hearing instrument is not in its T or M/T position. 2. Loop current (GAIN) has been set too low.	1. Change the hearing instrument to its T or M/T position. 2. Adjust it as required (see "Adjustment of loop current (GAIN)").
No sound.	1. Interruption in loop or poor contact at connection of loop to strip terminal (the green "Loop Current" lamp blinks or does not light up despite the signal reaching the TS-150). 2. No power (the red "Power" lamp does not light up).	1. Correct the interruption. 2. Check that the cord is correctly connected.
AGC lamp (yellow) lights up continuously or not at all.	Input level is too low or too high.	Adjust it as required (see "Adjustment of input levels").
Poor sound from the microphone	PHANTOM voltage is incorrectly installed.	Switch phantom voltage on/off (see "Microphone for the TS-150").
TS-150 overheats (its normal temperature is up to 60°C).	Product defect.	Pull out plug and contact the supplier.
The sound precedes the image.	The TV set's image processing makes the sound appear to precede the image.	Adjust the DELAY knob until the sound is synchronized.

AGC (Automatic gain control)

The TS-150 is equipped with an automatic system that controls the input signal. When the level control to the input is correctly adjusted, the AGC lamp should light up during the loudest parts of the talk or music. The TS-150 allows the time it takes for the signal to return to its normal level to be changed. When it has been correctly installed, strong signals are attenuated.

Supplier:

.....

Installation instructions for the TS-150


Preparation

- **Identify possible sources of interference:** switch on a hearing instrument, an inductive loop receiver or a field-strength meter and localise these sources. Any equipment causing interfering buzzing in the hearing instrument should be switched off or moved away where possible.
- **Signal sources:** typical sources are sound signals from TV sets, radios, music equipment and microphones. Only one signal source at a time should be switched on in addition to the microphone.
- **Placement:** the TS-150 should be mounted firmly on a wall. Avoid placement closer than 1 m to TV sets and potential interference sources. The TS-150 should be mounted in an area with good air circulation and should not be exposed to excessive heat from cookers and other electrical appliances.
- **The loop cable** should be run horizontally. It should be attached or glued firmly along the skirting board or edge of the ceiling (see “Notes”, point 1) around the room to be covered by the inductive loop system. The loop is laid in a pattern and with a cable type that assures matched loading (see below).
- **Cable type and area:** the loop cable may be of the single or multi-wire type. The thickness of the cable is selected for a total resistance of between 0.5 – 1 Ohm including the connecting cable (See “Selection of loop cable”).

Apart from the usual installation tools, signal cables and loop cables, the installer should be equipped with a field-strength meter and a telemagnetic receiver with headphones (e.g. Companion).

Mounting

1. The loop cable is attached or glued along the skirting board or ceiling edge around the area (room) to be covered. The loop cable is then connected to LOOP OUTPUT (see next page).
2. The amplifier is mounted on a wall with the mounting material supplied (drilling template, see centrefold).
3. The signal cables from the sound sources can be either plugged directly into the phono or Scart contacts or connected to the screw contacts on the connection block.


Connections


- RADIO-L: Signal input (e.g. from left channel of a radio or other sound system).
- RADIO-R: Signal input (e.g. from right channel of a radio or other sound system).
- VIDEO: Sound and video signal from/to video or DVD recorder (connected to the TV).
- TV: Sound and video signal from/to TV (connected to video or DVD recorder).
- LOOP OUTPUT: Pins 1 and 3 are outputs for the inductive loop. (Pin 2 is merely a connection point)
- Connection block: Pins 4 – 10 are alternative connection points for TV, RADIO-L & R, BAL.

Adjustments

After mounting, the signal levels must be adjusted.
(Check that the DIL switch for AGC = ON and SLOW and that PHANTOM = ON).

1. Switch off all sound sources.
2. Connect the TS-150 with the supplied equipment cable to a wall socket (220 V 50 Hz) and check that the red "POWER" LED lights up.
The yellow "AGC" LED usually also goes on, but goes off after 15-20 seconds.
3. Wait about 30 seconds until only the red "POWER" LED is on.


Adjusting the input levels

• RADIO inputs:

The trim potentiometer is shared by both radio inputs (L & R).

1. Adjust the trim potentiometer to its minimum (anticlockwise).
2. Switch on the sound source connected to RADIO-L and RADIO-R. Wait until the AGC lamp goes out if it went on when the sound source was switched on. Select a programme or station with a fairly constant level of speech or music and set the source loudspeaker to normal listening level.
3. Adjust RADIO slowly clockwise until the AGC lamp just goes on.
4. Switch off the sound source.

- **Balanced signals:**

If your sound source has a balanced signal, it should be connected to pins 8, 9 and 10 on the connection block: screen (ground) on pin 9 and the two signal leads on pins 8 and 10. The RADIO trim potentiometer often needs to be turned to its end stop (clockwise).

- **TV and VIDEO inputs:**

The TS-150 has a common trim potentiometer for TV and video inputs. Sound from the TV program and/or video recorder can reach the inductive loop by connecting signal cables from the TV set/video recorder to the TS-150 in two different ways:

- a) With a Scart cable from a Scart contact.
- b) With an RCA cable from the audio output.

a) From a Scart contact (EURO contact): the TS-150 is supplied with a Scart cable, about 2 meters long. A direct connection is made between the Scart contact (sometimes called the EURO contact) on the TV and the TV contact on the TS-150. If the Scart contact on the TV is connected to a video recorder, the video cable should be moved from the TV set to the Scart contact for VIDEO on the TS-150. An internal automatic program selector ensures that the sound and picture on the TV set correspond. The trim potentiometer for the TV on the TS-150 is set to its maximum when supplied. If the sound on the inductive loop is distorted or if the AGC lamp lights up continuously, it may be necessary to adjust the trim potentiometer slowly anticlockwise until the AGC lamp goes out and then clockwise until the AGC lamp just goes on.

b) From the audio output: As a rule, TV sets connected to satellite receivers, decoders, video recorders etc., have no free Scart contact, but are often equipped with an A/V output (three RCA contacts marked “Audio-L”, “Audio-R” and “Video”). This output is well suited as a sound source for the loop amplifier.

The TV audio outputs (often colour-coded red and white respectively) are connected to the TV inputs on the TS-150 **either** via a SCART–RCA/VHS INPUT adapter (see diagram on right), **or** both signals to pin 4 with the screen (ground) to pin 5 on the connection block.

SCART-RCA / VHS adapter

- i) Adjust the TV trim potentiometer on the TS-150 to its minimum (anticlockwise).
- ii) Connect the audio cables and switch the TV on.
- iii) If the AGC lamp goes on, you must wait until it goes off and then adjust the trim potentiometer slowly clockwise until the AGC lamp just goes on.

• About synchronization of sound (DELAY)


This should normally be in OFF mode. If the sound should need to be delayed, follow the instructions on page 3, "Synchronization of sound - The DELAY knob".

Microphone for the TS-150


Several different types of microphones may be used for the TS-150: dynamic, electret or capacitor microphones. The microphone input has a standard ¼-inch (2-pin) mono jack. A 10 V voltage (phantom measurement) is obtained on the contact if the DIL switch PHANTOM is ON. This voltage is necessary with the use of electret and capacitor microphones. Their properties and sensitivity can vary depending on their type and make. So it may be necessary to adjust the microphone input. The trim potentiometer for this purpose is marked MIC.

When using a:

- a. **Dynamic microphone:** Connect the microphone to the MIC input. Turn the DIL switch PHANTOM OFF. Switch on the microphone and speak into it with a normal voice. The microphone should be held at a distance of 5-10 cm from your mouth. The orange lamp MIC.PRIORITY should go on and remain on when you speak. If not, adjust it with the trim potentiometer MIC or speak more loudly.

- b. **Electret/capacitor microphone:** Adjust the trim potentiometer MIC down. Turn the DIL switch PHANTOM ON. Connect the microphone and switch it on. Speak into the microphone at a suitable distance depending on its make. Adjust MIC slowly clockwise until the orange MIC.PRIORITY lamp goes on.

Adjusting the sound

The TS-150 is set up for a flat frequency curve when supplied. Normally it is not necessary to change the base or treble setting on the loop amplifier, as this is adjusted in the hearing instrument. If required, however, these may be adjusted with the trim potentiometer TONE.

AGC or not?

The amplifier has an automatic gain control (AGC). It should ensure that the signals from the sound sources are kept within an acceptable level and prevent the signals on the loop becoming too strong or too weak. If the signals from the sound sources come from extra loudspeaker outputs, it is important to activate this function. If the signals to the TS-150 come from Scart, audio or line outputs, the AGC function is not necessary and can be switched off. Most hearing instruments have the AGC function connected in constantly to prevent the sound from the loop system becoming too loud for the ear.

Adjusting the loop current (GAIN)

To ensure good sound quality in the hearing instrument, it is important for the field strength to be correct and to comply with the requirements of the international standard IEC60118-4 (see Note 2 on the next page). The field strength depends on factors such as the loop current and must be adapted to the size of the room and its construction material. So it is important for the installation to be carried out correctly with the aid of a field-strength meter.

Switch on one of your sources and select a program with voice or music.

- a. Adjust the GAIN clockwise until the green LOOP CURRENT lamp lights up.
- b. Place the field-strength meter centrally in the loop (room) about 1.2 m above the floor and adjust the GAIN until the deflection of the field-strength meter oscillates around 0 dB. The oscillation should follow the signal from the sound source. The loudest parts of the speech/music should produce a deflection of up to +12 dB on the meter.
- c. Listen to the sound quality via the earphones of the field-strength meter or via a loop receiver (e.g. Companion).
- d. The user should now check the settings and evaluate the sound. The hearing instrument is placed in the T-position. During this evaluation, the user should sit on a chair in the centre of the room. *When the system is correctly installed, the sound level in the hearing instrument should be about equally strong from the inductive loop as from a person speaking at a distance of about 80 cm. It should not be necessary to adjust the sound level on the hearing instrument when changing its position from M to T.*

- e. When the system has been installed as stipulated by IEC60118-4, its quality should be checked for compliance with the values noted in the certificate (see “Certificate for inductive loop”). *For best results, the telecoil of the hearing instrument should be positioned vertically in the loop field. This is usually the case when the user sits and looks straight ahead. If he leans his head forwards or backwards, the sound usually becomes weaker. It also becomes weaker if he stands or sits directly above or in the immediate vicinity of the loop wire.*

Selection of loop cable

As already mentioned in the “Preparation” section, it is desirable for the resistance in the loop wire to be between 0.5 and 1 Ohm. This resistance becomes increasingly critical as the loop area gets larger. When the loop cable is selected according to the table below, the TS-150 normally satisfies the requirement in the IEC60118-4 standard as long as the rest of the settings are correct.

Wire area (mm ²), copper wire	Loop area (m ²)
0.25	Under 10
0.5	10 – 50
1 or 0.25 *	50 – 150

*For loop area of 50-150 m², a thinner cable can be used in order to decrease the temperature of the loop amplifier

Notes

- Ceiling-mounted loop:** If the loop is run along the ceiling, the ceiling height must not exceed 2.4 – 2.6 m. If it is greater, the loop must be mounted on the wall about 2.5 m above the floor.
- Extract from IEC60118-4:** Manufacturers of hearing instruments and loop amplifiers must observe the requirement of the standard in the development and production of the equipment. The loop system must be installed in such a way as to ensure maximum compliance with the standard. The standard recommends the following field strengths, measured at 1000 Hz in the centre of the loop (room) 1.2 m above the floor:
 - Nominal field strength : 0 dB (flat scale) = 100 mA/m.
 - Maximum field strength : 12 dB (flat scale) = 400 mA/m. at signal peaks (<0,125 s).
 - Frequency range : 100 Hz – 5000 Hz +/- 3 dB relative 1000 Hz.
 - Coverage area : More than 80% of loop area.
 - Signal/noise ratio : Less than -20 dB (A-weighted).

Technical data

Operating voltage	230V AC 50-60 Hz
Loop area	To IEC60118-4 up to 150 m ²
Frequency range	To IEC60118-4 at 100 – 5000 Hz ± 3 dB
Frequency range (+/-3dB)	90-15000Hz (delay off), 90-8500Hz (delay on)
Distortion	< 1%
Delay	Adjustable: off, 30-500 ms
DIL switch	
AGC ON - OFF	Connection/disconnection of AGC (automatic gain control)
FAST	Recovery time approx. 5 seconds
SLOW	Recovery time approx. 30 seconds
PHANTOM ON	+10 V to tip of microphone jack
PHANTOM OFF	No supply voltage to the microphone

Inputs	Clamp pin no.	Sensitivity AGC Off
TV	4 – 5 (5 is ground)	65 mV
RADIO-L	6 – 7 (7 is ground)	65 mV
RADIO-R	8 – 9 (9 is ground)	65 mV
BAL (Balanced)	8 – 10 (9 is ground)	40 mV

Via contacts	Contact type	
TV	Scart	70 mV
VIDEO	Scart	70 mV
RADIO-L	Phono (RCA)	65 mV
RADIO-R	Phono (RCA)	65 mV
MICROPHONE	¼-inch mono-jack	3 mV
	Preferred	Attenuates other signals

Loop connection		
LOOP OUTPUT	Clamp pins 1 – 3	Resistance 0.5 – 4 Ohm

Indicators (LED)	
POWER (Red)	Lights up continuously when TS-150 is connected to the power supply.
LOOP CURRENT (Green)	Lamp lights when current flows in the loop
AGC (Yellow)	Lamp should light when the input signal is in the AGC range
MIC. PRIORITY (Orange)	Lamp should light up when the microphone has voice input

Any issues relating to the EU Medical Device Directive 93/42/EEC should be directed to AB Transistor Sweden.

Cleaning, safety and maintenance

Cleaning and maintenance

Use a vacuum cleaner with a soft brush to clean the product on the outside. Use a slightly damp cloth to wipe it on the outside. Ensure that the cloth is not wet so that moisture cannot penetrate the amplifier. Never use cleaning products or solvents. From time to time check that the cables and connections are intact.

Repairs

If your product malfunctions, it must be repaired by a qualified technician. Do not attempt to open the case of the device since this would invalidate the warranty. If your product requires service, please contact your hearing care professional for assistance.


Waste disposal


This product contains electrical or electronic equipment and should be disposed of carefully in the interest of your safety and the environment. Please contact your local hearing care professional concerning disposal of the product.

General precautions

- Low-energy radiation (e.g. from burglar alarms, room surveillance systems, radio equipment, mobile telephones) will not damage the device. They could, however, momentarily affect the sound quality or create noise in the device.
- Warning: External devices connected to the product must be safe according the requirements of EN 60601-1-1, EN 60335-1 or EN 60065, as appropriate.

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MANUFACTURER:

in the sense of EU Medical Device Directive
93/42/EEC:

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Service and product-related questions: see
contact details for each country or
"International inquiries".