PolyMem®

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
PolyMem helps reduce spread of inflammatory reaction into surrounding, uninjured areas

(A) Incision only

(B) Incision with gauze

(C) Incision with PolyMem

This series of images shows the width of the spread of the inflammatory cells, in muscle, around an incision. The dark portion of the scale in each image (each segment is 100µm) represents the spread of the zone of the inflammatory reaction around the center line of the incision. In images A and B, there is no difference in the spread of the inflammatory reaction around the center of the injury. In image C, notice how PolyMem reduces the spread of the inflammation into the surrounding tissue. Statistically, PolyMem reduces the spread of the inflammatory reaction into the surrounding undamaged tissue by approximately 25 percent⁽⁸⁾.

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
Referenser:

1. Rodeheaver GT, Smith SL, Thacker JG, Edgerton MT, Edlich RF. Mechanical cleansing of contaminated wounds with a surfactant. *Am JSurg* 1975; 129(3):241-5.
2. Hokett SD, Cuenin MF, O'Neal RB et al. Pluronic polyol effects on human gingival fibroblast attachment and growth. *J Periodontol* 2000; 71(5):803-9.
3. Fluhr JW, Gloor M, Lehmann L, Lazzarini S, Distanti F, Berardesca E. Glycerol accelerates recovery of barrier function in vivo. *Acta Derm Venereol* 1999; 79(6):418-21.
4. Vandeputte J. Clinical Relevant Discoveries beyond occlusion when using a glycerin hydrogel. Poster Presented at 1995 Advanced Wound Care Symposium.
5. Lee BY, Carr RD, Herz B, Savino JA, Madden RA, Yadegar J. Polymeric Membrane Wound Care Dressing vs. Gauze (Standard Treatment). Scientific Exhibit. American College of Surgeons Annual Conference. October 20-32, 2003. Chicago, IL, USA.
6. Cutting KC, Weigand C, Vowden P. Wound inflammation and the role of dressings. *Wounds International* 2015; 6(2): 44-46.
7. European Wound Management Association. Position document: Pain and trauma at wound dressing changes. MEP Ltd: London, 2002.
8. Beitz, AJ, Newman A, Kahn AR, Ruggles T, Eikmeier L. A Polymeric Membrane Dressing with Antinociceptive Properties: Analysis with a Rodent Model of Stab Wound Secondary Hyperalgesia. *The Journal of Pain*. Feb 2004;5(1):38-47.
9. Knight KL. Chapter 3. Inflammation and Wound Repair In Cryotherapy in Sport Injury Management. Human Kinetics. 1995. Champaign, IL.
10. Merrick MA. Secondary injury after musculoskeletal trauma: a review and update. *Journal of Athletic Training* 2002;37(2): 209-217.
11. Clay CS, Chen WYJ. Wound pain: the need for a more understanding approach. *Journal of Wound Care*. April 2005;14(4):181-184.
12. Abraham SE. Pain Management in wound care. *Podiatry Management*. June/July 2006:165-168.
13. Wulf H, Baron R. The Theory of Pain in European Wound Management Association Position Document Pain at Wound Dressing Changes, Medical Education Partnership, London UK, 2002; page 8-11.
14. Levine JD, Reichling DB. Chapter 2 Peripheral Mechanisms of Inflammatory Pain. In Wall PD, Melzak R, Editors. *Textbook of Pain*. 4th edition. Edinburgh, UK: Churchill Livingstone, 1999. pages 59-84.
15. Fields HL. Chapter 1 Introduction & Chapter 2 The Peripheral Pain Sensory System In *Pain* New York: McGraw-Hill 1987 pages 1-40.
16. Brook I. Clinical clues to diagnosis of anaerobic infections <http://www.medscape.com/viewarticle/495997> accessed July 18, 2010.
17. Wolcott RD, Rhoads DD, Dowd SE. Biofilms and chronic wound inflammation. *Journal of Wound Care* 2008; 17(8):333-341.
18. Dow G, Browne A, Sibbald R.G. Infection in chronic wounds: Controversies in diagnosis and treatment. *Ostomy/Wound Management* 1999;45(8):23-40.
19. Foresman PA, Ethridge CA, Rodeheaver G. A wound healing evaluation on partial-thickness rat wounds. Symposium on Advances in Skin and Wound Care. 1991 Annual Meeting. Poster Presentation. Health Management Publication.

Individual product literature for specific indications and instructions for use. See:

PolyMem, PolyMem Silver, PolyMem Wic, PolyMem Wic Silver, PolyMem Max, PolyMem Max Silver, Shapes, Shapes by PolyMem, Ferris and design are trademarks of Ferris. The marks may be registered or pending in the US Patent and Trademark Office and in other countries. © 2018 Ferris Mfg. Corp. All rights reserved.