



If foam rubber is present at the attachment area we recommend that 4 cm. be removed so that **Flexiwise**® grips firmly onto the metal frame.

Om mjukgummi finns vid fästpunkten så rekommenderas att cirka 4 cm. av mjukgummit tas bort så att **Flexiwise**® får ett stadigt grepp om röret.



www.flexiwise.com

E-mail: info@flexiwise.com

Phone: +46 707 77 85 74

Copyright © 2010 Flexiwise AB, All rights reserved. Flexiwise ® is a national and international registered trademark of Flexiwise AB.