

VenoTrain[®] ulcertec

QUALITIES

The combination of an overstocking with a patented rhomboid knit and a skin-friendly liner for fixing the wound dressing in place makes the item easy to put on without any help from others.

The liners are also designed to provide a continuous low level of compression during the night, which can be individually ordered as a set of three and are suitable for washing at high temperatures.

- **easy to put on**
- **constant pressure**
- **discreet appearance**
- **X5+ sizing system**

COLORS

NATURAL

(Overstocking)

WHITE

(Liner)

MATERIAL CONTENT

39 – moderate 87 % polyamide
13 % elastane

46 – strong 69 % polyamide
31 % elastane

CARE INSTRUCTIONS FOR OVERSTOCKING



CARE INSTRUCTIONS FOR LINER





CEAP-
CLASSIFICATION

C0

C1

C2

C3

C4

C5

C6

LYMPHEDEMA
LIPOEDEMA

PREVENTION

STYLES


VenoTrain ulcertec OVERVIEW

| | |
|--------------------------|--|
| SIZING SYSTEM | X5+ (XS – XL) |
| COMPRESSION CLASS | 39 – moderate (34 – 46 mmHg) Liner with gentle compression 46 – strong (34 – 46 mmHg) Liner with moderate compression |
| TOE | Overstocking – open Liner – closed |

STANDARD SOLUTIONS

| | |
|------------------------|--|
| EDGING/TOP BAND |  AD |
| Top band | ● |

CUSTOM-MADE SOLUTIONS

| | |
|------------------------|--|
| EDGING/TOP BAND |  AD |
| Top band | ● |

VenoTrain ulcertec, 39 – moderate

- compression strength 34 – 46 mmHg
- Liner with gentle compression
- venous ulcer of the leg with good tendency to heal
- light liner for patients who can only tolerate a little compression at night

VenoTrain ulcertec, 46 – strong

- compression strength 34 – 46 mmHg
- Liner with moderate compression
- stagnating ulcers with a poor tendency to heal
- Liner with moderate compression providing wearing tolerance and a good acceptance level of compression at night