

# Advanced, premium-quality speech aid for a more natural expression

An electrolarynx is often used right after surgery, before healing enough to speak with a voice prosthesis. Many also choose to use an electrolarynx as a backup, or a temporary means of communication.

- To communicate immediately after surgery while healing
- To communicate while learning to speak with a voice prosthesis
- As a backup in case other form of speaking fails

## Provox® TruTone® EMOTE®

### Simple and easy to use

- Intuitive adjustments, no tools needed
- Easy-set modes for push-button simplicity
- Quick and easy monotone modes

Allows changes in voice volume from whispering to loud

Allows for variance in the tone of voice for a more natural expression

### Up to one month battery life and USB charging

- Charges just like your phone
- Most users experience 25-35 days between charges

Good speech intelligibility

Comfortable and ergonomic hourglass design



Product name	Qty.	Ref.no
Provox TruTone EMOTE	1 pc	7439

Includes: Extra-sharp sound-head cap, rechargeable battery, oral adapter, oral tubes, USB cable

### Contact us for more information

#### Head office:

Atos Medical AB,  
Hyllie Boulevard 17,  
SE-215 32 Malmö, Sweden

Tel: +46 (0) 415 198 00,  
Email: [info@atosmedical.com](mailto:info@atosmedical.com)

#### Manufacturer:

Atos Medical AB,  
Kraftgatan 8,  
SE-242 35 Hörby, Sweden