

# VenoTrain® ulcertec

## FOR VENOUS LEG ULCERS

The combination of an overstocking with a patented rhomboid knitted fabric and a skin-friendly liner for holding the wound dressing in place makes this compression stocking easy to put on without any need for assistance.

The liners are also ideal for providing a continuous low level of compression during the night. These are available to order separately in a pack of three and are suitable for washing at high temperatures.

### FEATURES

- Coordinated compression of overstocking and liner
- Unique compression profile that accelerates the healing process while ensuring exceptional patient comfort\*
- Differentiated treatment based on two designs within one compression class

#### 39 – moderate

- Compression class 3 (34–46 mmHg)
- Liner with gentle compression providing a compression pressure in the lower range of compression class 3 in combination with the overstocking

#### 46 – strong

- Compression class 3 (34–46 mmHg)
- Liner with moderate compression providing a compression pressure in the upper range of compression class 3 in combination with the overstocking

\* Jünger, M., Wollina, U., Kohnen, R., & Rabe, E. (2004). Efficacy and tolerability of an ulcer compression stocking for therapy of chronic venous ulcer compared with a below-knee compression bandage: results from a prospective, randomized, multicentre trial. *Current medical research and opinion*, 20(10), 1613–1623.

\*\* The AG liner additionally has a silicone top band on the upper edge of the stocking. It consists of polyamide, elastane, and silicone.

### COLORS

BEIGE  
OVERSTOCKING



WHITE  
LINER



**NEW**

Liner (moderate)  
now available as  
AG version

### DETAILS

**SIZING SYSTEM** X5+ (XS–XL)

#### COMPRESSION CLASSES

- 39–moderate: 34–46 mmHg  
Liner with gentle compression
- 46–strong: 34–46 mmHg  
Liner with moderate compression

#### TOE

- open (overstocking)
- closed (liner)

#### MATERIAL CONTENT

LINER AG\*\*

- 39–moderate: 68 % polyamide, 32 % elastane

LINER AD

- 39–moderate: 87 % polyamide, 13 % elastane
- 46–strong: 69 % polyamide, 31 % elastane

OVERSTOCKING

- 39–moderate / 46–strong: 79 % polyamide, 21 % elastane

#### CARE INSTRUCTIONS

- Overstocking:
- Liner:

#### ● STANDARD SIZES

#### ○ MADE-TO-MEASURE SOLUTIONS

#### TOP BAND EDGING



AD



AG

Top band



Silicone sensitive



C0

C1

C2

C3

C4

C5

C6