

Bauerfeind Nordic AB
 Avser positionsnr: 76, 77, 78, 79, 80, 81, 82, 83, 84, 85, 86, 87, 118, 119,
 120, 121, 122, 123, 124, 125, 126, 127, 128, 129

VenoTrain® micro

QUALITIES

With a particularly high microfiber content that makes it comfortable to wear, the VenoTrain micro is one of the most popular compression stockings. Available in many designs and a wide selection of colors, it highlights the individual style of your customers. Our individually colored batik versions in marine or black are particular eye-catchers.

- **50 % microfibers: breathable, kind to the skin and easy to clean**
- **figure-enhancing panty area (pantyhose) and comfort sole**
- **opaque full-fashioned stocking quality**
- **also available as a Design Edition with the patterns "Tango" and "Jive"**

MATERIAL CONTENT

- Ccl 1** 65 % polyamide
 35 % elastane
- Ccl 2** 55 % polyamide
 45 % elastane

CARE INSTRUCTIONS



COLORS AND DESIGN EDITION

	Tango	Jive	Batik Black	Batik Navy
CREAM				
HONEY				
CARAMEL				
AMARETTO				
ALMOND				
ESPRESSO				
BLACK				
ANTHRACITE				
NAVY				

VenoTrain-STYLIST

Select your favorite design quickly and easily:
www.bauerfeind.com/venotrain-stylist

As an added extra, all styles and colors of the VenoTrain micro can be optionally dyed over the top in navy or black using a batik technique. Delivered within 10 working days.

Colors may differ due to printing.





CEAP-
CLASSIFICATION

C0

C1

C2

C3

C4

C5

C6

LYMPHEDEMA
LIPOEDEMA





PREVENTION

STYLES





VenoTrain micro OVERVIEW

SIZING SYSTEM	PERFECT-FIT (S-XL)
COMPRESSION CLASS	Ccl 1 (18-21 mmHg) Ccl 2 (23-32 mmHg)
TOE	open, closed

STANDARD SOLUTIONS

EDGING TOP BAND	 AD	 AG	 AT	 ATU**
Top band	●		●	●
Lace sensitive*		●		
Silicone dots***	●	●		
Silicone dot pattern		●		

CUSTOM-MADE SOLUTIONS

EDGING TOP BAND	 AD	 AG	 AT	 ATU**
Top band	●		●	●
Lace sensitive*		●		
Silicone dots***	●	●		
Silicone dot pattern		●		

* From May 1, 2017, the following will apply to VenoTrain micro: The Lace Sensitive silicone top band will replace the Lace silicone top band. The Lace Sensitive silicone top band will be available only for the Design Edition and fashion colors until April 30, 2017.

** not available as Design Edition

*** not available as Design Edition, not in fashion colors